

[KO] Cytokeratin 7 Mouse mAb

ABZ1126

Product Overview

Name	[KO] Cytokeratin 7 Mouse mAb
Accession(Primary)	P08729
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human KRT7 protein fragments expressed in E.coli.
Calculated MW	Calculated MW: 51 kDa; Observed MW: 55 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P08729)

Gene Name	KRT7
About Gene Symbol	KRT7 (Keratin 7) is a keratin cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.