

## [KO] Cytokeratin 8 Mouse mAb

ABZ1127

### Product Overview

|                             |                                                                           |
|-----------------------------|---------------------------------------------------------------------------|
| <b>Name</b>                 | [KO] Cytokeratin 8 Mouse mAb                                              |
| <b>Accession(Primary)</b>   | P05787                                                                    |
| <b>Host Species</b>         | Mouse                                                                     |
| <b>Isotype</b>              | IgG2a                                                                     |
| <b>Immunogen</b>            | Purified recombinant human KRT8 protein<br>fragments expressed in E.coli. |
| <b>Species Reactivity</b>   | Rat                                                                       |
| <b>Calculated MW</b>        | Calculated MW: 54 kDa; Observed MW: 54 kDa                                |
| <b>Recommended Dilution</b> | WB: 1/500-1/1000; IHC: 1/50-1/100; ICC:<br>1/100-1/200                    |

### Target information (P05787)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | KRT8 |
| <b>About Gene Symbol</b> |      |