

[KO] Cytokeratin 18 Mouse mAb

ABZ1128

Product Overview

Name	[KO] Cytokeratin 18 Mouse mAb
Accession(Primary)	P05783
Host Species	Mouse
Isotype	IgG2b
Immunogen	Purified recombinant human KRT18 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P05783)

Gene Name	KRT18
About Gene Symbol	