

[KO] L1CAM Mouse mAb

ABZ1129

Product Overview

Name	[KO] L1CAM Mouse mAb
Accession(Primary)	P32004
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human L1CAM produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 140 kDa; Observed MW: 250 kDa
Recommended Dilution	WB: 1/500-1/1000; FC: 1/50-1/100

Target information (P32004)

Gene Name	L1CAM
About Gene Symbol	L1CAM (L1 Cell Adhesion Molecule) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.