

[KO] Galectin 3 Mouse mAb

ABZ1131

Product Overview

Name	[KO] Galectin 3 Mouse mAb
Accession(Primary)	P17931
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human LGALS3 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/100-1/200

Target information (P17931)

Gene Name	LGALS3
About Gene Symbol	LGALS3 (Galectin 3) is an extracellular matrix laminin protein involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.