

[KO] Leukotriene A4 Hydrolase Mouse mAb

ABZ1135

Product Overview

Name	[KO] Leukotriene A4 Hydrolase Mouse mAb
Accession(Primary)	P09960
Host Species	Mouse
Isotype	IgG2b
Immunogen	Full length human recombinant protein of human LTA4H produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 69 kDa; Observed MW: 69 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P09960)

Gene Name	LTA4H
About Gene Symbol	LTA4H (Leukotriene A4 Hydrolase) is a protease enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.