

[KO] HIF1 beta Mouse mAb

ABZ1136

Product Overview

Name	[KO] HIF1 beta Mouse mAb
Accession(Primary)	P27540
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human ARNT produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 87 kDa; Observed MW: 110 kDa
Recommended Dilution	WB: 1/500-1/1000; ICC: 1/100-1/200; FC: 1/50-1/100

Target information (P27540)

Gene Name	ARNT
About Gene Symbol	ARNT (Aryl Hydrocarbon Receptor Nuclear Translocator) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.