

[KO] c-Met Mouse mAb

ABZ1138

Product Overview

Name	[KO] c-Met Mouse mAb
Accession(Primary)	P08581
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human MET protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 156 kDa; Observed MW: 156 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/100-1/200

Target information (P08581)

Gene Name	MET
About Gene Symbol	