

[KO] Ki67 Rat mAb

ABZ1139

Product Overview

Name	[KO] Ki67 Rat mAb
Accession(Primary)	P46013
Host Species	Rat
Isotype	IgG1
Immunogen	Purified recombinant human MKI67 protein fragments expressed in E.coli.
Calculated MW	Calculated MW: 359 kDa; Observed MW: 359 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P46013)

Gene Name	MKI67
About Gene Symbol	