

[KO] Nibrin Mouse mAb

ABZ1140

Product Overview

Name	[KO] Nibrin Mouse mAb
Accession(Primary)	O60934
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human NBN protein fragments expressed in E.coli.
Calculated MW	Calculated MW: 85 kDa; Observed MW: 105 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (O60934)

Gene Name	NBN
About Gene Symbol	