

[KO] NM23A Mouse mAb

ABZ1143

Product Overview

Name	[KO] NM23A Mouse mAb
Accession(Primary)	P15531
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human NME1 produced in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P15531)

Gene Name	NME1
About Gene Symbol	NME1 (NME/NM23 Nucleoside Diphosphate Kinase 1) is a protein kinase and signaling regulator involved in cell migration and invasion. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.