

[KO] NNMT Mouse mAb

ABZ1145

Product Overview

Name	[KO] NNMT Mouse mAb
Accession(Primary)	P40261
Host Species	Mouse
Isotype	IgG2b
Immunogen	Full length human recombinant protein of human NNMT produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; FC: 1/50-1/100

Target information (P40261)

Gene Name	NNMT
About Gene Symbol	