

[KO] NONO Mouse mAb

ABZ1146

Product Overview

Name	[KO] NONO Mouse mAb
Accession(Primary)	Q15233
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human NONO protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 65 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q15233)

Gene Name	NONO
About Gene Symbol	