

[KO] EXOSC1 Mouse mAb

ABZ1147

Product Overview

Name	[KO] EXOSC1 Mouse mAb
Accession(Primary)	Q9Y3B2
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human EXOSC1 produced in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q9Y3B2)

Gene Name	EXOSC1
About Gene Symbol	EXOSC1 (Exosome Component 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.