

[KO] MRPS7 Mouse mAb

ABZ1148

Product Overview

Name	[KO] MRPS7 Mouse mAb
Accession(Primary)	Q9Y2R9
Host Species	Mouse
Isotype	IgG2a
Immunogen	Purified recombinant human MRPS7 protein fragments expressed in E.coli.
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/100-1/200

Target information (Q9Y2R9)

Gene Name	MRPS7
About Gene Symbol	