

[KO] PFKP Mouse mAb

ABZ1149

Product Overview

Name	[KO] PFKP Mouse mAb
Accession(Primary)	Q01813
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human PFKP produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 86 kDa; Observed MW: 86 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; FC: 1/50-1/100

Target information (Q01813)

Gene Name	PFKP
About Gene Symbol	PFKP (Phosphofructokinase, Platelet) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.