

[KO] PHKG2 Mouse mAb

ABZ1150

Product Overview

Name	[KO] PHKG2 Mouse mAb
Accession(Primary)	P15735
Host Species	Mouse
Isotype	IgG2a
Immunogen	Purified recombinant human PHKG2 protein fragments expressed in E.coli.
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P15735)

Gene Name	PHKG2
About Gene Symbol	