

[KO] PIK3C2A Mouse mAb

ABZ1151

Product Overview

Name	[KO] PIK3C2A Mouse mAb
Accession(Primary)	O00443
Host Species	Mouse
Isotype	IgG2a
Immunogen	Purified recombinant human PIK3C2A protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 191 kDa; Observed MW: 191 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (O00443)

Gene Name	PIK3C2A
About Gene Symbol	PIK3C2A (Phosphatidylinositol-4-Phosphate 3-Kinase Catalytic Subunit Type 2 Alpha) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic targ