

## [KO] PMS2 Mouse mAb

ABZ1152

### Product Overview

|                             |                                                                            |
|-----------------------------|----------------------------------------------------------------------------|
| <b>Name</b>                 | [KO] PMS2 Mouse mAb                                                        |
| <b>Accession(Primary)</b>   | P54278                                                                     |
| <b>Host Species</b>         | Mouse                                                                      |
| <b>Isotype</b>              | IgG2b                                                                      |
| <b>Immunogen</b>            | Full length human recombinant protein of human PMS2 produced in 293T cell. |
| <b>Calculated MW</b>        | Calculated MW: 96 kDa; Observed MW: 150 kDa                                |
| <b>Recommended Dilution</b> | WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/100-1/200                        |

### Target information (P54278)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | PMS2 |
| <b>About Gene Symbol</b> |      |