

[KO] PUS7 Mouse mAb

ABZ1153

Product Overview

Name	[KO] PUS7 Mouse mAb
Accession(Primary)	Q96PZ0
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human PUS7 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	WB: 1/500-1/1000 FC: 1/50-1/100
Recommended Dilution	WB: 1/500-1/1000; FC: 1/50-1/100

Target information (Q96PZ0)

Gene Name	PUS7
About Gene Symbol	