

[KO] CHCHD3 Mouse mAb

ABZ1154

Product Overview

Name	[KO] CHCHD3 Mouse mAb
Accession(Primary)	Q9NX63
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human CHCHD3 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q9NX63)

Gene Name	CHCHD3
About Gene Symbol	CHCHD3 (Coiled-Coil-Helix-Coiled-Coil-Helix Domain Containing 3) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.