

[KO] PPP2R1B Mouse mAb

ABZ1155

Product Overview

Name	[KO] PPP2R1B Mouse mAb
Accession(Primary)	P30154
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human PPP2R1B protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 66 kDa; Observed MW: 66 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P30154)

Gene Name	PPP2R1B
About Gene Symbol	