

[KO] PPP5C Mouse mAb

ABZ1156

Product Overview

Name	[KO] PPP5C Mouse mAb
Accession(Primary)	P53041
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human PPP5C produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P53041)

Gene Name	PPP5C
About Gene Symbol	PPP5C (Protein Phosphatase 5 Catalytic Subunit) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.