

[KO] OTUB1 Mouse mAb

ABZ1157

Product Overview

Name	[KO] OTUB1 Mouse mAb
Accession(Primary)	Q96FW1
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human OTUB1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 31 kDa; Observed MW: 36 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q96FW1)

Gene Name	OTUB1
About Gene Symbol	OTUB1 (OTU Deubiquitinase, Ubiquitin Aldehyde Binding 1) is a deubiquitinase involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellul