

[KO] PKA RII alpha Mouse mAb

ABZ1158

Product Overview

Name	[KO] PKA RII alpha Mouse mAb
Accession(Primary)	P13861
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human PRKAR2A produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 55 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P13861)

Gene Name	PRKAR2A
About Gene Symbol	PRKAR2A (Protein Kinase CAMP-Dependent Type II Regulatory Subunit Alpha) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research