

[KO] MAPK1 Mouse mAb

ABZ1159

Product Overview

Name	[KO] MAPK1 Mouse mAb
Accession(Primary)	P28482
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human MAPK1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P28482)

Gene Name	MAPK1
About Gene Symbol	MAPK1 (Mitogen-Activated Protein Kinase 1) is a receptor tyrosine kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.