

[KO] ERK5 Mouse mAb

ABZ1160

Product Overview

Name	[KO] ERK5 Mouse mAb
Accession(Primary)	Q13164
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human MAPK7 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 88 kDa; Observed MW: 110 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (Q13164)

Gene Name	MAPK7
About Gene Symbol	