

[KO] JNK2 Mouse mAb

ABZ1161

Product Overview

Name	[KO] JNK2 Mouse mAb
Accession(Primary)	P45984
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human MAPK9 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P45984)

Gene Name	MAPK9
About Gene Symbol	