

[KO] MEK1 Mouse mAb

ABZ1162

Product Overview

Name	[KO] MEK1 Mouse mAb
Accession(Primary)	Q02750
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human MAP2K1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q02750)

Gene Name	MAP2K1
About Gene Symbol	MAP2K1 (Mitogen-Activated Protein Kinase Kinase 1) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.