

[KO] MEK2 Mouse mAb

ABZ1163

Product Overview

Name	[KO] MEK2 Mouse mAb
Accession(Primary)	P36507
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human MAP2K2 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P36507)

Gene Name	MAP2K2
About Gene Symbol	MAP2K2 (Mitogen-Activated Protein Kinase Kinase 2) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.