

[KO] LXN Mouse mAb

ABZ1164

Product Overview

Name	[KO] LXN Mouse mAb
Accession(Primary)	Q9BS40
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human LXN produced in 293T cell.
Calculated MW	Calculated MW: 26 kDa; Observed MW: 35 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q9BS40)

Gene Name	LXN
About Gene Symbol	LXN (Latexin) is a protease enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.