

[KO] OTUD7B Mouse mAb

ABZ1166

Product Overview

Name	[KO] OTUD7B Mouse mAb
Accession(Primary)	Q6GQQ9
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human OTUD7B protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 93 kDa; Observed MW: 93 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (Q6GQQ9)

Gene Name	OTUD7B
About Gene Symbol	