

[KO] Twinfilin 1 Mouse mAb

ABZ1167

Product Overview

Name	[KO] Twinfilin 1 Mouse mAb
Accession(Primary)	Q12792
Host Species	Mouse
Isotype	IgG2b
Immunogen	Purified recombinant human TWF1 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q12792)

Gene Name	TWF1
About Gene Symbol	