

[KO] PTP1B Mouse mAb

ABZ1168

Product Overview

Name	[KO] PTP1B Mouse mAb
Accession(Primary)	P18031
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human PTPN1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P18031)

Gene Name	PTPN1
About Gene Symbol	