

[KO] SHP2 Mouse mAb

ABZ1169

Product Overview

Name	[KO] SHP2 Mouse mAb
Accession(Primary)	Q06124
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human PTPN11 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Recommended Dilution	WB: 1/500-1/1000; FC: 1/50-1/100

Target information (Q06124)

Gene Name	PTPN11
About Gene Symbol	PTPN11 (Protein Tyrosine Phosphatase Non-Receptor Type 11) is a protein phosphatase involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.