

[KO] Bax Mouse mAb

ABZ1170

Product Overview

Name	[KO] Bax Mouse mAb
Accession(Primary)	Q07812
Host Species	Mouse
Isotype	IgG2b
Immunogen	Purified recombinant human BAX protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q07812)

Gene Name	BAX
About Gene Symbol	BAX (BCL2 Associated X, Apoptosis Regulator) is a protein kinase and signaling regulator involved in apoptosis. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.