

[KO] PEX5 Mouse mAb

ABZ1171

Product Overview

Name	[KO] PEX5 Mouse mAb
Accession(Primary)	P50542
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human PEX5 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 71 kDa; Observed MW: 80 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (P50542)

Gene Name	PEX5
About Gene Symbol	