

[KO] Rad23B Mouse mAb

ABZ1172

Product Overview

Name	[KO] Rad23B Mouse mAb
Accession(Primary)	P54727
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human RAD23B protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 50 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P54727)

Gene Name	RAD23B
About Gene Symbol	RAD23B (RAD23 Homolog B, Nucleotide Excision Repair Protein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.