

[KO] Rb Mouse mAb

ABZ1173

Product Overview

Name	[KO] Rb Mouse mAb
Accession(Primary)	P06400
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human RB1 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 106 kDa; Observed MW: 106 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P06400)

Gene Name	RB1
About Gene Symbol	RB1 (RB Transcriptional Corepressor 1) is a protein phosphatase involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.