

[KO] S100 alpha 6 Mouse mAb

ABZ1176

Product Overview

Name	[KO] S100 alpha 6 Mouse mAb
Accession(Primary)	P06703
Host Species	Mouse
Isotype	IgG2b
Immunogen	Full length human recombinant protein of human S100A6 produced in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 10 kDa; Observed MW: 10 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (P06703)

Gene Name	S100A6
About Gene Symbol	