

[KO] S100P Mouse mAb

ABZ1177

Product Overview

Name	[KO] S100P Mouse mAb
Accession(Primary)	P25815
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human S100P produced in 293T cell.
Calculated MW	Calculated MW: 10 kDa; Observed MW: 10 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/100-1/200

Target information (P25815)

Gene Name	S100P
About Gene Symbol	