

[KO] ATG3 Mouse mAb

ABZ1178

Product Overview

Name	[KO] ATG3 Mouse mAb
Accession(Primary)	Q9NT62
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human ATG3 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 42 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q9NT62)

Gene Name	ATG3
About Gene Symbol	