

[KO] SH3GL1 Mouse mAb

ABZ1179

Product Overview

Name	[KO] SH3GL1 Mouse mAb
Accession(Primary)	Q99961
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human SH3GL1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 50 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; FC: 1/50-1/100

Target information (Q99961)

Gene Name	SH3GL1
About Gene Symbol	SH3GL1 (SH3 Domain Containing GRB2 Like 1, Endophilin A2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.