

[KO] Fascin Mouse mAb

ABZ1180

Product Overview

Name	[KO] Fascin Mouse mAb
Accession(Primary)	Q16658
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human FSCN1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q16658)

Gene Name	FSCN1
About Gene Symbol	