

[KO] SOD2 Mouse mAb

ABZ1181

Product Overview

Name	[KO] SOD2 Mouse mAb
Accession(Primary)	P04179
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human SOD2 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 25 kDa; Observed MW: 20 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (P04179)

Gene Name	SOD2
About Gene Symbol	SOD2 (Superoxide Dismutase 2) is a metabolic enzyme involved in cellular metabolism. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.