

[KO] SPR Mouse mAb

ABZ1182

Product Overview

Name	[KO] SPR Mouse mAb
Accession(Primary)	P35270
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human SPR produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; FC: 1/50-1/100

Target information (P35270)

Gene Name	SPR
About Gene Symbol	