

[KO] BRAF Mouse mAb

ABZ1183

Product Overview

Name	[KO] BRAF Mouse mAb
Accession(Primary)	P15056
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human BRAF produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 105 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P15056)

Gene Name	BRAF
About Gene Symbol	