

[KO] TRAP alpha Mouse mAb

ABZ1184

Product Overview

Name	[KO] TRAP alpha Mouse mAb
Accession(Primary)	P43307
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human SSR1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 32 kDa; Observed MW: 37 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P43307)

Gene Name	TRAPA
About Gene Symbol	SSR1 (Signal Sequence Receptor Subunit 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.