

[KO] STAT1 Mouse mAb

ABZ1185

Product Overview

Name	[KO] STAT1 Mouse mAb
Accession(Primary)	P42224
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human STAT1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 87 kDa; Observed MW: 87 kDa
Recommended Dilution	WB: 1/500-1/1000; ICC: 1/100-1/200

Target information (P42224)

Gene Name	STAT1
About Gene Symbol	STAT1 (Signal Transducer And Activator Of Transcription 1) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cell