

[KO] STAT3 Mouse mAb

ABZ1186

Product Overview

Name	[KO] STAT3 Mouse mAb
Accession(Primary)	P40763
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human STAT3 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 88 kDa; Observed MW: 88 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P40763)

Gene Name	STAT3
About Gene Symbol	STAT3 (Signal Transducer And Activator Of Transcription 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.